

# ComfyCouch C102

For Aphaeresis and Blood Donation



Introducing the Comfy Couch C102 for Safety and Comfort of the Donor and Convenience of the Technician. Aesthetically designed with focus on quality, safety and hygiene. Intelligently designed and proven to reduce incidence of vasovagal attack associated with blood donation.

**QUALIMED**  
*Assuring Quality*



## Superior Design

Position can be shifted to shock position within 10 seconds of actuation. Single profile construction enables the donor profile to be maintained unchanged during conversion to shock position, ensuring minimal discomfort to the donor. The wedge shaped armrest promotes comfortable vena puncture facilitating optimal blood flow. A Large resting surface enables comfortable seating and long foot section for easy foot support for taller donors.



## Safety

The Unique position of the foot during donation ensures adequate blood supply to the head as per hemodynamics of the human body to minimize the incidence of "Vasovagal attack". Controlling buttons on the remote are logically placed and does not confuse. High stability ensures no tipping over up to 150kgs. Provision for all devices and accessories required for blood donation and height adjustment of the couch ensures safety from RSI [Repetitive Strain Injuries] for staff.



## Hygiene

Single profile construction ensures that all surfaces are smooth and easy to clean. The upholstery is not stitched and do not provide any opportunity for moisture or liquids to lodge. Smooth, dirt resistant design of the remote without any sharp edges or dirt traps in the construction. Accessories trays can be removed for easy cleaning. Easy to clean transparent vinyl cover that protects foot area.



➤ *Height can be adjusted by 200mm from minimum*

**Unique 4 way Action Hand Rest**

- Height Adjustable.
- Can be raised or swung out for easy entry of Donor.
- Angle of Arm rest can be adjusted for Optimum Venipuncture.
- Angle to profile can also be set.



Multiple recline position with quick conversion to trendelenburg position in 10 seconds.

Comfortable donor position with long foot profile provides comfort to even the tallest Donor. Special transparent vinyl covering at foot section avoids dirt accumulation and prevents damage by footwear.

Smooth Fiberglass housing and hence internal motor and hydraulics not exposed

## TECHNICAL DATA:

Input power	220 ± 10% ~ VAC, 50 Hz.
Power	120 W
Output Socket Voltage	220±10% V, 50Hz
Output Socket Power	60W Max
Maximum Load	150 kg.
Seat Width	520mm
Profile length	1920mm
Seat height	620mm
Width of arm rest	150mm
Length of arm rest	480mm
Wheel diameter	100mm

## COLOUR OPTIONS



Navy Blue



Steel Grey



Maroon



Beige

\*Any other colour as per customer requirement if order qty is more than 10units

Qualimed is a specialist manufacturer of Couches for Aphaeresis and Blood Donation.

**COMFY COUCH C101** is the base model ideally suited for a modern blood donation center. Comfy Couch C101 also comes with an option for fixing brackets instead of wheels for installation in Mobile Blood Bank Bus/Van

**COMFY COUCH D4** is an advanced model with 4 motors by which the seat, back rest, leg rest and height can be adjusted. This is ideally suited for long term procedures like aphaeresis, dialysis and oncology treatment.



COMFY COUCH C101



COMFY COUCH D4